

Differences Between a Wheelchair and a Transport Chair



Wheelchair

- Large 18"-24" rear wheels, Small front wheels. (38-88lbs)
- Rotating handrail on rear wheels allows user to push him/herself.
- Heavy duty construction for long duration sitting with padded armrests, backrest & seat.
- Push-lock mechanism feature locks wheels in place for safely entering and exiting chair.
- Accessories available - trays, cup-holders, bags, cushions, etc.



Transport Chair

- 8" or 12" rear wheels, small front wheels. (25-32lbs)
- Requires assistance of a second person to push from behind and to engage handle brakes when needed.
- Travel friendly and lightweight; compact design easily folds closed.
- Narrow frame to pass through most doors and aisle ways.
- Seatbelt attached for safety.

